

MAGELLAN^ε

NANOPARTICLE TRACE CHARACTERIZATION IN WATER

*The ultimate nanoparticle analyzer
for water characterization*



Push your analysis few steps further !

IDEAL FOR

Online monitoring of water treatment
Membrane integrity monitoring
Ultrapure water monitoring
Water quality insurance
Filtration processes
Pollution detection
... and more

www.cordouan-tech.com



- Unique patented technology
- Nanoparticle size and concentration analysis
- Trace analysis down to ppt (ng/L)

Breakthrough technology inside


MAGELLAN^ε is a unique **nanoparticle analyzer** based on a patented breakthrough technique called Laser Induced Breakdown Detection. **MAGELLAN^ε** measures **size distribution and concentration** of **nanoparticle traces** in water with an unprecedented sensitivity and resolution.

MAGELLAN^ε is the result of a successful joint development between **CORDOUAN Technologies** and the Karlsruhe Institute of Technology (KIT) in Germany.



LIBD Principle

LIBD is a field proven technology issued from more than 20 years of research in KIT-INE; basically, a nanosecond pulsed laser is focused into a liquid containing nanoparticles to be detected. Each time a particle crosses the laser beam, a **plasma** is created and detected by an **acoustic wave sensor**.

Size distribution and concentration are then deduced from the plasma statistic by the use of an advanced proprietary algorithm and calibration data.




[Download technical notes
www.cordouan-tech.com](http://www.cordouan-tech.com)



- Broadest particle size measurement
- Sensitivity: more than 10 000 times superior to conventional light scattering techniques


Unprecedented performances, Unequaled sensitivity & resolution



Nanoparticle trace characterization in water

Magellan's Key benefit

- ✓ Number density : as low as 10^4 particles/mL
- ✓ High size resolution through 8 channels
- ✓ Compact and robust design : transportable and operational in all environments
- ✓ Versatile and easy-to-use :
 - Variable flow rate
 - No sample preparation required
 - Limited maintenance (highly reliable design)
- ✓ Broadest particle size range : from 10 nm to $1\mu\text{m}$
- ✓ Unequaled sensitivity : down to ng/L (ppt) range
- In line and under pressure measurements (up to 60 bars)
- Advanced Proprietary algorithm for accurate results
- Various sample cell configurations: static or flow-through cell



Targeted applications

Drinking and industrial water quality insurance



Membrane integrity online monitoring (fouling & breakthrough)



Pollution detection of natural water resources




**Push your analysis
few steps further**

Colloids in primary coolant of nuclear reactor system characterization



Ultrapure water monitoring for semiconductor industry



MAGELLAN^ε

Nanoparticle trace characterization in water

Specifications

Particle size range	10 – 1000 nm (Up to 8 channels)
Concentration range	10 ⁴ to 10 ¹¹ part/ml (Typically 10 ppt with 20nm particles)
Accuracy	+/-10% (depending on measurement time)
Calibration	Certified reference materials polystyrene (NIST) or customer reference
Sample preparation	None – Standard laboratory practices for ultra-trace sample analysis
Standard sample cell configuration	Static cell : 3.5ml / Flow-through cell: up to 5 bars – 275µL - 4 mL/min (pump depending)
On demand sample cell configuration	In situ / High pressure cell: up to 60 bars - 4 mL/min

Signal processing

Measurement technology	Laser Induced Breakdown Detection (LIBD)
Laser source	DPSS - Power <50mW - Pulsed laser 100 Hz – TEM00 - @532 nm – Lifetime 2 Giga shots
Data processing algorithm	Proprietary algorithm software developed with KIT : ASTROLIBD [®]
Detector	Piezoelectric acoustic wave sensor (Patented)

General


Operating conditions / Storage conditions	5°C to 40°C / -10°C to 50°C – Relative humidity < 70% non condensing
Computer interface	USB 2.0 – RJ45 – Serial Port - Windows 7 professional – 64-bit
Dimensions	40 cm x 30 cm x 40 cm (HWD)
Weight	Analyzer < 25 kg / Computer < 10kg
Power	100-240 VAC, 50-60 Hz, 10A max

System Compliance

CE certification	CE marked product - Class I laser product, EN 60825-1:2001, CDRH
Ingress Protection (IP)	Housing IP54, Optics IP65 – (protection from dust and water)



Flow-through cell



Standard static cell

The **MAGELLAN** analyzer and its components are protected by:
German Patents No.: DE 19833339C1
 DE 19602048C2, DE 102006051227B9;
European Patent No.: EP 1 918 694 A1;
Worldwide Patent No.: WO 0006993A1;
 other patent pending.



Contact

sales@cordouan-tech.com
 service@cordouan-tech.com
 11, avenue de Canteranne
 33600 Pessac – France
 Phone : +33 (0)556 158 045
 Fax : +33 (0)547 747 491

www.cordouan-tech.com

 @CordouanTech